



PerkinElmer[™]
life sciences.

Sample Preparation & Liquid Scintillation Counting

*From Professionals
For Professionals*

J.Thomson
Worldwide C&S Product Manager

What is an LSC Cocktail ?

An LSC Cocktail can be considered as :-

Solvent (Aromatic)

Scintillators (Fluorescing Agents)

This mixture is sometimes called a

_Lipophilic Cocktail



Solvent

Scintillators

Surfactant (Detergent)

This mixture is also called an

Emulsifying Cocktail.




Why is an LSC Cocktail used ?

The majority of radioactive species are present in an aqueous form, and as such are not miscible with aromatic solvents.

The presence of surfactants (detergents) in the cocktail enables an aqueous sample to come into intimate contact with the aromatic solvent by forming a stable microemulsion.

How does an LSC cocktail work ?



How does an LSC cocktail work ?

Scintillation Counting

1

Liquid Scintillation Cocktail



Components:

Solvent: Typically toluene, xylene pseudocumene, or an alkyl benzene type solvent.

Emulsifier: A detergent type molecule that ensures proper mixing of aqueous samples.

Fluor: A fluorescent solute.

Process: Radioactive Sample is added to scintillation cocktail.




How does an LSC cocktail work ?




Scintillation Counting

2

Beta particles are emitted, which cause solvent molecules to become excited.



The energy of the solvent molecules is transferred to the fluor molecules, which in turn emit light.



Why are these components in an LSC Cocktail?

The aromatic solvent is necessary since it contains a high density of π electrons, necessary for the efficient transfer of the energy of radioactive decay.

The surfactant is needed to enable a stable microemulsion to be formed when aqueous samples are present, necessary for stable conditions over the counting period.

The scintillators are present to emit a light pulse which is within the optimum detection range of the photomultiplier tube.

What Aromatic Solvents are in LSC Cocktails ?

1. Toluene



CAS No. 108-88-3

RPH = 100




2. Xylene (mixed isomers)



CAS No. 1330-20-7

RPH = 110




3. Pseudocumene (1,2,4-trimethylbenzene)



CAS No. 95-63-6

RPH = 112




4. Dodecyl Benzene (LAB -Linear Alkyl Benzene)



CAS No. 123-01-3

RPH = 94




5. Di-isopropylnaphthalene (DIN)



CAS No. 38640-62-9

RPH = 112




6. 1-Phenyl-1-Xylyl Ethane (PXE)



CAS No. 6196-95-8

RPH = 110



What Detergents are used in LSC Cocktails ?


1. Ethoxylated Alkylphenols
(Non-Ionic detergent)
CAS No. 9016-45-9



2. Mono-/Di- phosphate ester
(Anionic detergent)
CAS No. 298-07-7



3. Sodium di-octylsulphosuccinate
(Anionic detergent)
CAS No. 577-11-7




What Scintillators are used in LSC Cocktails ?

Primary Scintillators

A. 2,5-diphenyloxazole (**PPO**)
(Primary Scintillator)

$C_{15}H_{11}NO$


CAS No. 92-71-7



B. **Butyl PBD** [2-(4-tert-butylphenyl)-5-(4-biphenyl)-1,3,4-oxadiazole]
(Primary Scintillator)

$C_{24}H_{22}N_2O$

CAS No. 15082-28-7



What Scintillators are used in LSC Cocktails ?


Secondary Scintillators

C. 1,4-Bis(2-methylstyryl)benzene (**bis-MSB**)

(Secondary Scintillator)

$C_{24}H_{22}$

CAS No. 13280-61-0




D. 1,4-Bis(4-methyl-5phenyl-2oxazolyl)benzene (**Dimethyl POPOP**)

(Secondary Scintillator)

$C_{26}H_{20}N_2O_2$

CAS No. 3073-87-8



What Vials are best?

??
?
?
?
?
?
?
?
?



?
?
?
?
?
?
?
?
?

Glass is used because :

- ✓ It provides unparalleled optical clarity (good visibility)
- ✓ It is chemically inert, making it suitable for use with aggressive reagents and when solvents are used.

Plastic is used because :

- ✓ It exhibits lower background levels than glass.
- ✓ It is combustible and therefore easier for waste disposal.
- ✓ It is shatterproof and therefore safer in the laboratory.

What size of vials are available ?

Vials (glass and plastic) are available in sizes ranging from **<4mL up to 20mL.**


This maximum of 20mL is fixed due to the dimensions of currently used photomultiplier tubes (2 inch diameter).

Why is a sample presented in a liquid form ?

Liquid samples are preferred due to the fact that the scintillation counting process takes place most efficiently in an LSC cocktail (due to 4π geometry).


Since certain liquids are miscible and can therefore mix together rapidly and completely, samples in a liquid medium are preferred.

4π Geometry (Homogeneous)



$$\text{Surface Area} = 4\pi r^2$$

2π Geometry (Heterogeneous)



$$\text{Surface Area} = 2\pi r^2$$

Inorganic Anions

Chlorides, Sulfates, Phosphates, Nitrates, etc.

Potential Obstacles

Di- and Tri- valent anions sometimes cause instability.
Both volume and concentration dependent.
Some metallic salts form colored solutions.

Remedy

- ✓ Select a cocktail suitable for use with these difficult anions (**Ultima Gold and HiSafe series**)
- ✓ Dilute sample with water or increase cocktail volume.
- ✓ Use a cocktail which is resistant to color quenching (**Ultima Gold and HiSafe series**)

Mineral Acids and Organic Acids

Hydrochloric Acid, Nitric Acid, Sulfuric Acid, Phosphoric Acid, TCA, Acetic Acid etc.

Potential Obstacles

Certain strong mineral acids interact with the cocktail components causing color development and performance changes. Conc. sulfuric acid will cause color formation and sulfonation. TCA can result in acid-induced luminescence.

Remedy

- ✓ Do not add concentrated acids to the cocktail.
- ✓ Always dilute the acid with water.
- ✓ Use a cocktail which has proven performance with acidic samples (**Ultima Gold AB & LLT**).

Inorganic Alkalis

Sodium Hydroxide, Potassium Hydroxide, etc.

Potential Obstacles

Alkalis usually cause chemiluminescence.
Prolonged contact can result in color formation.
Conc. potassium hydroxide will give high "backgrounds".

Remedy

- ✓ Use a cocktail with proven resistance to CLM
(**Hionic-Fluor**)
- ✓ Dilute the alkali to reduce the effect.
- ✓ Acidify the cocktail / sample mixture.
- ✓ Do not store for prolonged periods.
- ✓ Change to Sodium Hydroxide; KOH contains ^{40}K .

Potential Obstacles:

Quench effect of organic solvent i.e.

Nitro > Sulphide > Halide > Amine > Ketone >
Aldehyde > Organic Acid > Ester > (Water) >
Alcohol > Ether > Hydrocarbon

Remedy:

1. Change the organic solvent to a less quenching alternative and/or use a more chemical quench resistant cocktail.
2. Always check to ensure that the cocktail accepts the mixture.

1. **Direct Counting** (addition of raw sample to a suitable cocktail) is suitable for :
Urine, Serum, Plasma and Water soluble protein
2. **Solubilization** (dissolution/digestion of sample) is suitable for :
Blood, Plasma, Serum, Feces, Homogenates & Bacteria
3. **Combustion** (oxidation/burning sample to convert the organic carbon and hydrogen to $^{14}\text{CO}_2$ and $^3\text{H}_2\text{O}$ respectively) is suitable for all of the above sample types without exception.

Urine :-

Color quench:

Reduce sample volume.

Use correct cocktail.

Precipitation of proteins:

Arrange to count samples within 24 hours. Correct cocktail selection (**Pico-Fluor MI, Ultima Gold LLT**)

Additional information

Urine from small animals more difficult

Recommended LSC cocktails for Urine

Cocktail	Type	Max. urine* capacity @20°C	Recommended loading
Ultima Gold LLT	Safer	10.0mL	2.0mL
Ultima Gold	Safer	8.0mL	2.0mL
Pico-Fluor MI	Classical	4.0mL	4.0mL

* Results based on human urine

Note :-

1. It is possible to add up to 1mL of either ethanol or isopropyl alcohol per 10mL cocktail, to suppress the appearance of the protein precipitate.
2. Urine samples from small animals are usually difficult. With these, either keep sample volume low, or dilute the urine with water.

Other LSC cocktails for Urine

Cocktail	Type	Max. capacity at 20°C	Comments
Ultima Gold	Safer	8.0mL	Slight wispy precipitate appears after 24hours at >2.0mL sample.
Ultima Gold XR	Safer	10.0mL	Slight wispy precipitate appears after 24hours at >3.0mL sample.
Ultima Gold AB	Safer	10.0mL	Slight wispy precipitate appears after 24hours at >2.0mL sample.
Ultima Gold LLT	Safer	10.0mL	Slight wispy precipitate appears after 24hours at >2.0mL sample.
Opti-Fluor	Safer	9.0mL	Slight wispy precipitate appears after 24hours at >3.0mL sample.
Emulsifier Safe	Safer	2.5mL	Slight wispy precipitate appears after 24hours at >3.0mL sample.
Pico-Fluor MI	Classical	4.0mL	Remains clear up to 4.0mL sample
Hionic-Fluor	Classical	1.5mL	>2.0mL sample forms a gel
Insta-Gel Plus	Classical	1.0mL	
Pico-Fluor 40	Classical	3.0mL	Slight wispy precipitate appears after 24hours at >3.0mL sample.

Problems with Biological Samples

Plasma/Serum :-

Sample volume:

Keep sample volume to <1mL per 10mL cocktail

Color quench:

Reduce sample volume.

Use correct cocktail (**Pico-Fluor MI, Ultima-Flo M**)

Precipitation of proteins:

Arrange to count samples within 24 hours.


Correct cocktail selection.

(Try adding 10% ethanol or IPA to inhibit protein precipitation) Plasma/Serum from small animals is more difficult.

Biological Sample



Solubilizer



Alkaline Systems

- ✓ Quaternary ammonium hydroxide, NaOH
- ✓ Method of action is by ***alkaline hydrolysis***.

Acidic Systems


- ✓ Perchloric Acid, Nitric Acid, Perchloric/Nitric Acid Mixture, Perchloric Acid/Hydrogen Peroxide Mixture
- ✓ Method of action is by ***acidic oxidation***.

Other Systems


Sodium Hypochlorite

Method of action is by ***oxidative bleaching***.


What happens during Solubilization?




1. **Add sample** to glass vial.




2. **Add solubilizer** and **heat** at $\sim 55^{\circ}\text{C}$.



3. **Add peroxide** to decolorize.



4. **Add LSC Cocktail**.



5. **Dark adapt** before **counting**.

Solubilizers & Recommended Cocktails



Soluene-350

- ✓ Organic-based solubilizer, formulated in toluene
- ✓ Classical industry-standard solubilizer
- ✓ For biological and plant samples
- ✓ Corrosive & Flammable

Recommended cocktail - **Hionic-Fluor**




SOLVABLE

- ✓ **Aqueous-based solubilizer (safer)**
- ✓ Can substitute for Soluene-350 in most cases
- ✓ For biological and plant samples
- ✓ Corrosive, non-flammable

Recommended cocktail - **Ultima Gold**

Good sample preparation
=
Good LSC counting

Health & Safety



The correct term to use with LSC cocktails is **SAFER**, i.e. safer than the tradition or classical solvents such as Toluene, Xylene and Pseudocumene.

NO chemical can be properly classified as SAFE.

- **Classical Cocktails** - Toluene, Xylene, Pseudocumene
 - ✓ Hazardous by inhalation
 - ✓ Hazardous by skin absorption
 - ✓ Irritating to skin and eyes
 - ✓ High vapor pressure
 - ✓ Low flash point (flammable)

- **Newer Safer Cocktails** - DIN, PXE, LAB
 - ✓ Low toxicity
 - ✓ Classified as harmless (for transport & storage)
 - ✓ Low vapor pressure
 - ✓ High flash point (non-flammable)
 - ✓ Biodegradable

Safety Phrases (Mandatory)



***Clearly visible
hazard warnings
& safety phrases***

PerkinElmer Safer Cocktails:

Testing for biodegradability is according to Test Method OECD 301E (equivalent to ISO Method 7827 1984). To pass, the sample must show >70% DOC removal within 28 days.

A cocktail which passes this test can be described as “**Readily Biodegradable**”.

Disposal Considerations

Waste products may be disposed of by controlled incineration. National and local legislation concerning the disposal of waste materials must be complied with.

“We do not advocate drain disposal of cocktail waste. Biodegradability information on Packard cocktails is provided to ensure that in the event of a small spillage, any waste that ends up in the sewerage system will be biodegraded.”

A few final thoughts !!!!!

***Any LSC cocktail will work with any LSC -
The trick is to get the best combination.***

Poor sampling can only produce a bad sample.

A badly prepared sample can only produce bad results.

***No amount of instrument sophistication can ever
produce good results from a badly prepared sample.***