

Analysis of ^3H & ^{14}C using the Raddec Pyrolyser

Ian Croudace & Phil Warwick



Geosciences
Advisory
Unit

RADDEC^{LTD}
Radiochemical and Decommissioning Solutions Ltd

Introduction

- The Pyrolyser series[©] of furnaces allow the simultaneous decomposition of up to 6 samples and the efficient oxidation of the liberated ^3H and ^{14}C species
- The furnace design has been thoroughly tested and proven over many years.
- Continued R&D ensures that upgrades become incorporated into the system design

Features


- 4 or 6 samples decomposed simultaneously in ~4 hours using a programmable thermal ramp.
- 10g Pt-alumina catalyst per tube allowing 20+ determinations.
- Rapid cool-down of sample zone enabling good cycle times between runs.
- Bubbler-traps (20mls) have >95% efficiency.
- The Pyrolysers have been in regular use for extraction of ^3H and ^{14}C from a wide range of materials (vegetation, fish, soil, sediment, concrete, metal etc).
- The design and proven effectiveness follows several years of testing of samples from intercomparison exercises, environmental studies and nuclear site decommissioning programmes.

Number of H-3/C-14 analysis

- 2002 H-3 (242)/C-14 (79)
- 2003 H-3 (376)/C-14 (127)
- 2004 H-3 (847)/C-14 (565)
- 2005 H-3 (996)/C-14 (153)

Max. number of samples that can be run using a single 6-tube furnace is 1352



(52 x 24 = 1352 p.a.)



Geosciences
Advisory
Unit

RADDEC^{LTD}
Radiochemical and Decommissioning Solutions Ltd

Purpose-built tube furnace for H-3/C-14 analysis (Pyrolyser 6™)



Geosciences
Advisory
Unit

RADDEC^{LTD}
Radiochemical and Decommissioning Solutions Ltd


Sample types

- Environmental samples

- Soil/Sediment, Fruit, Water, Grass, Milk, Fish, Sludge etc.

- Decommissioning samples

- Concrete, Brick, Asbestos, MMMF, Metal, Plastic, Desiccants, Paper, Electric wire, Sewage sludge, Graphite, Paint etc.




Geosciences
Advisory
Unit

RADDEC^{LTD}
Radiochemical and Decommissioning Solutions Ltd

Sample sizes

- Samples with low organic contents ;
Soil/Sediment, concrete, brick, metal etc. (1 - 10g)
- Samples with high organic contents ;
Biota, plastics, organic rich sediment/soil etc.
(Normally 0.5 - 1g but can be increased with a slow ramp)



Geosciences
Advisory
Unit

RADDEC^{LTD}
Radiochemical and Decommissioning Solutions Ltd

Counting efficiency and furnace recovery for H-3/C-14

- Technical information


	SQPE*	Counting efficiency	Furnace recovery**	Limit of detection***
H-3	719 - 729	18 – 20 %	> 90 %	0.020 Bq/g
C-14	720 - 750	66 – 73 %	> 95 %	0.015 Bq/g

* Typical SQPE of concrete samples

** Average value of 70 measurements using an organic H-3/C-14 standard

*** Using 5g sample size and 2 hour counting time on Quantulus

$$L_D(Bq/g) = \frac{2.71 + 4.65\sqrt{C}}{t} \times \frac{100}{E} \times \frac{100}{R} \times \frac{1}{M_g}$$




Geosciences
Advisory
Unit

RADDEC^{LTD}
Radiochemical and Decommissioning Solutions Ltd

Intercomparison results*

Sample	H-3/C-14 type	Measured value Bq/g $\pm 1\sigma$	Reference value Bq/g $\pm 1\sigma$
NPL (2002)	Tritiated water	20 \pm 10	20.04 \pm 0.18
NPL (2004)	Tritiated water	0.536 \pm 0.042	0.539 \pm 0.006
NPL (2004)	C-14 carbonate	23.4 \pm 1.6	24.4 \pm 0.5
FSA	H-3 thymidine (Milk)	4.72 \pm 0.66	Mean = 4.04 0.18 – 4.93
FSA	H-3 thymidine (Plaice)	4.42 \pm 0.30	Mean = 4.67 2.7 – 8.3

* There are no commercially available certified reference material for H-3/C-14 activities



Geosciences
Advisory
Unit

RADDEC^{LTD}
Radiochemical and Decommissioning Solutions Ltd

LSC⁺ -Plus



Enhanced LSC
data processing software

..... better and faster

Results for the H-3 in water (8ml + 12ml Gold Star™)

Report date : 7-Nov-2003

Customer : Raddec Ltd

Job reference : Raddec 372

Date samples received : 6-Nov-2003

Date of analysis : 7-Nov-2003

Working instruction number: - Raddec/RC/2022

Calibration report number: - Raddec/CAL/16


Counter S/N	Laboratory S/N	Reference date	H-3 Bq/ml	2 s.d.
1 INST STD		7-Nov-2003	23.224	2.548
2 6-169-1	372-1	7-Nov-2003	< 0.006	
3 6-169-2	372-2	7-Nov-2003	< 0.007	
4 6-169-3	372-3	7-Nov-2003	< 0.006	
5 6-169-4 STD		7-Nov-2003	0.115	0.014

All results are in Bq/ml and are decay corrected to the reference date (12.50 Years half life)


"< values" are limits of detection as defined by Currie, 1968

Uncertainties are at the 2 s.d. confidence level and are based on propagated method uncertainties

Analyst : A Other



- Calculates activities directly from counter files
- Eliminate transcription errors
- Calculates LODs (Currie)
- Built in Quality Control
- Statistical analysis of results
- Full diagnostic report
- Range of input file formats for all LSC counters




Geosciences
Advisory
Unit

RADDEC^{LTD}
Radiochemical and Decommissioning Solutions Ltd

Organisations currently using the Raddec Pyrolyser

- GAU (Geosciences Advisory Unit): 3 furnaces
 - AMEC-NNC (NIRAS) Ltd.: 5 furnaces
 - Scientifics Ltd.: 1 furnace
 - HPA (Health Protection Agency): 1 furnace (Formerly known as NRPB)
 - AWE plc. (Atomic Weapons Establishment): 1 furnace
 - Babcock BES : 1 furnace
 - DSTL (Defence Science & Technology Laboratory): 1 furnace
 - KAERI (South Korea) : 1 furnace
 - 4 installations pending
- 14 Pyrolyser furnaces are being used for both commercial and research purposes in the UK.
 - The Pyrolyser was invented, developed and commercialised by Raddec. The Pyrolyser design is copyright protected.



Geosciences
Advisory
Unit

RADDEC^{LTD}
Radiochemical and Decommissioning Solutions Ltd

Raddec Ltd
PO Box 1210
Southampton
SO17 1ZQ
U.K.

Tel: (+44) 07739 898344

Tel & Fax: (+44) 02380 559367

Email: sales@raddec.com

www.raddec.com