

Monte Carlo modeling of an electronic brachytherapy source using MCNP5 and EGSnrc




Stephen D. Davis and Larry A. DeWerd

University of Wisconsin, Madison, WI

13th UK Monte Carlo User Group meeting
March 29, 2007




Xoft Axxent™ electronic brachytherapy source



Photos courtesy of Xoft, Inc.



Xoft Axxent™ electronic brachytherapy source



X-Ray Probe Tip Detail

Photos courtesy of Xoft, Inc.



Motivation

- Brachytherapy uses air-kerma strength, S_K , as the measurement of source strength
- Air kerma rates can be measured using free-air ionization chambers, but S_K is defined *in vacuo*
- Correction is required for attenuation in air from the source to the point of measurement (100 cm in our case)
- Air cross sections change rapidly at photon energies < 20 keV
- Photon spectrum needs to be accurately known for a low uncertainty on the air attenuation correction




Photon spectra


- Measurements with high purity germanium spectrometer
 - Requires correction for the energy response of the detector
- Monte Carlo simulations
 - Assumes Monte Carlo code has the low energy physics included to model x-ray spectra accurately




Photon spectra at 178 cm in air



Photon spectra at 178 cm in air



Photon spectra at 0 cm in air



EGSnrc

- January 2005 release included model by Iwan Kawrakow for electron impact ionization for K- and L-shells with binding energies above 1 keV
- October 2005 release included C++ class library with general purpose geometry modeling and user code `cavity.cpp`
- February 2007 release allowed for scoring of photon spectra in a circular plane using `cavity.cpp`



EGSnrc vs. MCNP5

- Two codes have differences in:
 - Electron impact ionization
 - Binding effects and Doppler broadening for Compton interactions
 - Atomic relaxation
 - Electron transport
 - Variance reduction techniques
 - Scoring options




Doppler broadening


- “Broomstick problem” from MCNP5 LANL document
- 100 keV photons incident on 10 cm graphite rod with 10^{-6} cm radius



No Doppler broadening



With Doppler broadening




Atomic relaxation


- 30 keV photons incident on 10^{-5} cm thick yttrium slab




Atomic relaxation (cont.)




Atomic relaxation (cont.)



Atomic relaxation (cont.)



cavity.cpp model of Axxent™ source




Dose rates at 1 cm in water

MC code	$P(r_0, \theta_0)$ [cGy·h ⁻¹ ·μA ⁻¹]
MCNP5	345 ± 0.1%
EGSnrc	346 ± 2.0%
EGSnrc/MCNP5	1.003 ± 2.0%




Air kerma rates at 100 cm in air


MC code	$P(r_0, \theta_0)$ [cGy·h ⁻¹ ·μA ⁻¹]
MCNP5	$4.81 \cdot 10^{-2} \pm 0.1\%$
EGSnrc	$4.72 \cdot 10^{-2} \pm 3.1\%$
EGSnrc/MCNP5	$0.981 \pm 3.1\%$




Calculated photon spectra at 1 cm in vacuum




Calculated photon spectra at 1 cm in vacuum



Calculated photon spectra at 1 cm in vacuum



Calculated photon spectra at 1 cm in vacuum



Conclusions

- Both codes produce similar water dose rates, air kerma rates, and photon spectra, and the small differences can be explained by different treatment of the low energy physics
- Further work will be necessary to compare these results to measurements, and to establish the required air attenuation corrections with low uncertainties



Acknowledgements

- Xoft, Inc.
- Iwan Kawrakow
- NIST
- UWMRRC staff and students
- UW ADCL customers
- Grid Laboratory of Wisconsin (GLOW)

