


Digital Gamma-ray Spectroscopy & Imaging with Semiconductor detectors

Frontiers of gamma-ray spectroscopy





Dr Andy Boston
ajboston@liv.ac.uk




THE UNIVERSITY
of LIVERPOOL



Without Compton suppression shields



Compton continuum.
=> Large peak to total ratio





With BGO shielding



Less solid angle coverage
=> Big drop in efficiency


With highly segmented detectors





Path of γ -ray reconstructed to form full energy event
=> Compton continuum reduced
=> Excellent efficiency **~50%** @1MeV
=> Greatly improved **angular resolution** ($\sim 1^\circ$) to reduce **Doppler effects**

$$E'_\gamma = \frac{E_\gamma}{1 + (E_\gamma / mc^2)(1 - \cos \theta)}$$

AGATA : Two candidate configurations



Ge crystals size:
length 90 mm
diameter 80 mm




Digital Spectroscopy and Imaging

180 hexagonal crystals	3 shapes
60 triple-clusters	all equal
Inner radius (Ge)	22 cm
Amount of germanium	310 kg
Solid angle coverage	80 %
Singles rate	~50 kHz
6480 segments	
Efficiency: 40% ($M_\gamma=1$)	25% ($M_\gamma=30$)
Peak/Total: 65% ($M_\gamma=1$)	50% ($M_\gamma=30$)

120 hexagonal crystals	2 shapes
40 triple-clusters	2 shapes
Inner radius (Ge)	17 cm
Amount of germanium	220 kg
Solid angle coverage	74 %
Singles rate	~70 kHz
4320 segments	
Efficiency: 38% ($M_\gamma=1$)	21% ($M_\gamma=30$)
Peak/Total: 63% ($M_\gamma=1$)	47% ($M_\gamma=30$)

The "Standard" Germanium Shell

Idealized configuration to determine maximum attainable performance



$R_i = 15 \text{ cm}$
 $R_o = 24 \text{ cm}$
 230 kg of Ge

$$M_\gamma = 1 \rightarrow \epsilon_{\text{ph}} = 65\%$$


$$P/T = 85\%$$

$$M_\gamma = 30 \rightarrow \epsilon_{\text{ph}} = 36\%$$

$$P/T = 60\%$$

Assuming 5 mm Position Resolution

A high multiplicity event




27 gammas detected -- 23 in photopeak
 16 reconstructed -- 14 in photopeak



AGATA Design




3 different asymmetric hexagonal shapes are used





Triple cluster modular units in a single cryostat




The **AGATA demonstrator**: 5 triple clusters, 540 segments. Scheduled for completion 2008



2π of complete array




Completed array (6480 segments) with support structure




AGATA 1st symmetric cluster

Digital Spectroscopy and Imaging



Canisters AGATA and EUROBALL




AGATA

EUROBALL

Automated Scanning Tables

Liverpool System



- Parker linear positioning table
- Pacific scientific stepper motors
- 920MBq ^{137}Cs
- 1 collimator
- Singles/coincidence system



Precise position
calibration

“How well your basis fits your real data”


Liverpool Detector Scanning System




Data recorded using 14 bit 80Mhz digital electronics to disk.

Pulse shape analysis of detector response as a function of position.

Charge pulse response









Cs-137 Singles Fine scan









- 920MBq Cs-137 source
- 60 sec per position
- 1mm steps
- 1mm diameter collimator

- Max count rate = 720 cps
- Background = 40 cps
- 600keV CFD threshold
- Accepted triggers = 120 cps
- 180GB pre-sorted data


Core 662keV – Ring gated



Core 662keV & Fold 1, Core T30




Core 662keV & Fold 1, Core T90





Coincidence scan: set up

Physical segmentation
depth:






Collimation gap centred on
(w.r.t. crystal base):



Coincidence scan: set up




Which grid to choose ?

		advantages	drawbacks
	r, θ	cst. values of $t_{10-90\dots}$ cylindrical	not homogenous
	\sqrt{r}, θ	cst. values of $t_{10-90\dots}$ homogen., cylindr.	not the same x/y accuracy
	x, y, z	homogenous simple	not cylindrical large distances to grid
	hexagon, z	cylindrical compact	not compact in z not homogenous
	hexagonal compact	cylindrical maximum compactity	less "standard"
	Adaptive grid	optimum conditioning of the problem	not homogenous


Cartesian Grid

- Different colors show active regions for the different segments




New Quasi-cylindrical Grid

- Different colors show active regions for the different segments
- Spacing is equidistant in sensitivity




New Quasi-cylindrical Grid

- Different colors show active regions for the different segments
- Spacing is equidistant in sensitivity




Coincidence scan: plan





Segment pulses



Centre contact pulses





6.0
mm




15.5
mm




28.5
mm



45.5
mm




Electric Field Simulations : MGS




• Electric field simulations have been performed and details comparisons have been made with experimental pulse shape data.

agata_vertex geometry




Potential Mapping





The SmartPET Project




THE UNIVERSITY
of LIVERPOOL



Andy Boston
ajboston@liv.ac.uk



SmartPET December 2005



SmartPET detector depth response

“superpulse” pulse shapes for ^{137}Cs events versus depth






Image Charge Response

Image charge asymmetry varies as a function of lateral interaction position

- Calibration of asymmetry response

$$\text{Asymmetry} = \frac{\text{Area}_{\text{left}} - \text{Area}_{\text{right}}}{\text{Area}_{\text{left}} + \text{Area}_{\text{right}}}$$






Image Charge Analysis



One approach to applying simple PSA is to split this distribution equally into five regions

Image Charge Analysis


Asymmetry distribution as a function of scan (collimator) position (AC05)




Reconstructed Images

- Simple PSA techniques applied event-by-event
- Filtered Back Projection – ^{22}Na source

No PSA




PSA





Compton Camera

- $10\mu\text{Ci } ^{152}\text{Eu}$
- 60mm from SPET1
- Source rotated
- Zero degrees in 15° steps up to 60°
- Detector separation
- 3 - 11cm in 2cm steps
- Gates set on energies
- 779, 1408keV
- 2 ^{22}Na sources at different x and y





Analytical method



J Gillam, Monash



Iterative Method



J Gillam, Monash

Compton Camera measurements (Ge/Ge)


$E = 1408 \text{ keV}$



J Gillam, Monash

Digital Gamma-ray Spectroscopy & Imaging with Semiconductor detectors

Frontiers of gamma-ray spectroscopy



Dr Andy Boston
ajboston@liv.ac.uk



THE UNIVERSITY
of LIVERPOOL